



Discover
the Magic of
easyPhaco®

oertli®
S W I T Z E R L A N D



Fluidics on...

New and faster:
Oertli **easyPhaco**® –
the technology which makes
fluidics to your best friend

Turn up vacuum (500 mmHg/50 ml for peri-
staltic, 500 mmHg or more for venturi), and
let the elaborate fluidics concept of the
Oertli system work for you!

Even though it sounds implausible – Oertli
easyPhaco® Technology brings to you

- unprecedented chamber stability
- perfect emulsification
- efficient fragment aspiration

And all of this without the
undesired side effects hitherto
caused by high vacuum.




Oertli **easyPhaco**® –
the leading technology in
your hands

New and better:
The Oertli **easyTip**[®] 2.2 mm

Intelligent needle design and drastically improved fluidics properties – Oertli easyPhaco[®] technology brings visible and perceptible advantages.

- no turbulences and efficient aspiration without repulsion of fragments
- concentrated axial ultrasound energy delivery for perfect emulsification
- no vacuum surge thanks to a capillary aspiration path

The Oertli **easyTip**[®] 2.2 mm is designed for Oertli[®] OS3 and SwissTech systems with Hexadisq[®] hand piece. Together, the unparalleled system for more safety and better results.



It's better – it's different – it's Oertli

oertli®
S W I T Z E R L A N D

Oertli Instrumente AG
CH-9442 Berneck, Switzerland
Phone +41 (0)71 747 42 00
Fax +41 (0)71 747 42 90
www.oertli-instruments.com
info@oertli-instruments.com