


# Microsurgical


## DCR Sets



- OSD20A** DCR Set Angled 20g
- OSD23A** DCR Set Angled 23g
- OSDCRAP** DCR Set Angled Pediatric

Stainless steel probes. 4.5cm ~ 11cm length, with 45 degree bend 15 mm from tip. Solid silicone tubing 40 cm length. Silicone tubing goes inside through the steel probes which prevents disconnection during manipulations.

10 Per box




- OSD20S** DCR Set Straight 20g
- OSD21S** DCR Set Straight 21g
- OSD22S** DCR Set Straight 22g
- OSD23S** DCR Set Straight 23g
- OSD24S** DCR Set Straight 24g
- OSD25S** DCR Set Straight 25g
- OSD26S** DCR Set Straight 26g
- OSD27S** DCR Set Straight 27g
- OSD28S** DCR Set Straight 28g
- OSD29S** DCR Set Straight 29g
- OSD30S** DCR Set Straight 30g
- OSDCRSP** DCR Set Straight Pediatric

Stainless steel probes. 4.5cm ~ 11cm length, solid silicone tubing 40cm length. Silicone tubing goes inside through the steel probes which prevents disconnection during manipulations.

10 Per box

## Lacrimal Intubation Sets

21g 23g 25g 27g




- OSL21C** Lacrimal Intubation Set (Crawford) 21g
- OSL23C** Lacrimal Intubation Set (Crawford) 23g
- OSL25C** Lacrimal Intubation Set (Crawford) 25g
- OSL27C** Lacrimal Intubation Set (Crawford) 27g
- OSLESCP** Lacrimal Intubation Set (Crawford) Pediatric

Stainless steel probes with olive shape tips 11cm ~ 20cm length Silicone tubing - 35 cm length. Silicone tubing goes inside through the steel probes which prevents disconnection during manipulations.

5 Per box

21g 23g 25g 27g




- OSL21** Lacrimal Intubation Set 21g
- OSL23** Lacrimal Intubation Set 23g
- OSL25** Lacrimal Intubation Set 25g
- OSL27** Lacrimal Intubation Set 27g
- OSLESP** Lacrimal Intubation Set Pediatric

Stainless steel probes 11cm ~ 20cm length Silicone tubing - 35 cm length. Silicone tubing goes inside through the steel probes which prevents disconnection during manipulations.

5 Per box

20g 21g 23g 25g 27g



- OSLX20** Lacrimal Intubation Set (Crawford) 20g
- OSLX21** Lacrimal Intubation Set (Crawford) 21g
- OSLX23** Lacrimal Intubation Set (Crawford) 23g
- OSLX25** Lacrimal Intubation Set (Crawford) 25g
- OSLX27** Lacrimal Intubation Set (Crawford) 27g
- OSLXP** Lacrimal Intubation Set (Crawford) Pediatric

Stainless steel probes with olive shape tips and retrieval device. Probes 11cm ~ 20cm length Silicone tubing - 35 cm length. Silicone tubing goes inside through the steel probes which prevents disconnection during manipulations.

5 Per box

## Choppers

Phaco Chopper



**BM100** Phaco Chopper

1.4 mm tip with blade offset at 45 degree angle. Overall length 115 mm. Used to break up the nucleus.

5 Per box



Nucleus Chopper



**BM101** Nucleus Chopper

45 degrees, 6mm from tip. Tip angled 90 degrees, 1.5mm from end. Used as a second instrument to break up the nucleus within the capsular bag during two-handed phaco techniques.

5 Per box



## Manipulators

Buratto's Hook



**BM102** Buratto's Hook

Angled, 0.25mm diameter, blunt tip. Overall length 120mm. Used to manipulate the nucleus.

5 Per box



Micro Manipulator



**BM110** Micro Manipulator

Angled, 0.35mm diameter, blunt tip. Overall length 120mm. Used to manipulate the nucleus.

5 Per box



Lens Manipulator (Lester)



**BM111** Lens Manipulator (Lester)

Angled, 0.44mm diameter hourglass shape tip. Overall length 115mm. Used for IOL manipulation.

5 Per box



## Nucleus Manipulator



### BMI12 Nucleus Manipulator

Flattened paddle tip, 1.10mm tip. Overall length 115mm. Used to manipulate the nucleus.  
5 Per box



## IOL Manipulator (Jaffe-Maltzman)



### BMI20 IOL Manipulator (Jaffe-Maltzman)

Angled 45 degrees, 12mm from end, blunt tip. Fits .20mm diameter holes and larger.  
5 Per box



## IOL Manipulator (Sinskey)



### BMI21 IOL Manipulator (Sinskey)

Straight, blunt tip. Fits .20mm diameter holes and larger.  
5 Per box



## IOL Manipulator (Sinskey)



### BMI22 IOL Manipulator (Sinskey)

Angled 45 degrees, 15mm from end, blunt tip. Fits .20mm diameter holes and larger.  
5 Per box



## IOL Manipulator (Sinskey)



### BMI23 IOL Manipulator (Sinskey)

Blunt tip, formed 10mm from end. Fits .20mm diameter holes and larger.  
5 Per box



## Button Tip Manipulator



### BMI24 Button Tip Manipulator

Angled 45 degrees, .50mm diameter button tip. 15mm from end. May be used for IOL manipulation or iris retraction.

5 Per box



## Nucleus Rotator



### BMI113 Nucleus Rotator

Angled, 0.70mm diameter flat, blunt tip. Overall length 115mm. Used to manipulate the nucleus.

5 Per box



## Nucleus Rotator



### BMI125 Nucleus Rotator

45 degrees, blunt forked tip. 9mm from end. For rotating nucleus during phacoemulsification.

5 Per box



## LASIK Flap Roller



### BMI130 LASIK Flap Roller

14mm sleeve. Used after realignment. Acts as a squeegee to remove excess fluid from the interface, smooths flap and removes folds.

5 Per box



## LASIK Cap / Flap Elevator




### BMI131 LASIK Cap / Flap Elevator

Designed to automatically lift corneal cap/flap during corneal lamellar surgery. Angled end is used for separating epithelium around previous incision. Spatulated end used to lift corneal flap.


5 Per box




**Capsule Retractor**




**BM132 Capsule Retractor**  
 Angled, 0.60mm diameter semi-ball. Overall length 115mm. Used to retract capsular bag or manipulate nucleus.  
 5 Per box



**Micro Retractor**




**BM133 Micro Retractor**  
 30 degrees. Micro size allows iris or capsula retraction without interfering with other surgical procedures.  
 5 Per box




## Spatulas


**Spatula**




**BM150 Spatula**  
 .80mm (45 degrees). 10mm from end. Smooth round edges allow for safe manipulation within the eye.  
 5 Per box




**Spatula**





**BM151 Spatula**  
 .80mm. Smooth round edges allow for safe manipulation within the eye.  
 5 Per box



**Orucoglu Spatula**



**BM155 Orucoglu Spatula®**  
 LASIK spatula with two ends. One end with a blunt and curved design especially for easy entry and separation of flap. Modified Sinskey end for separation and elevation of the new flap.  
 5 Per box


Punctal Dilator and Probe (Khan)



**BM160** Punctal Dilator and Probe (Khan)

25mm medium taper, blunt. Wilder style dilator mated with a size 0 Bowman lacrimal probe. .70mm x 45mm. Dual function instrument allows for single handed transition from dilating to probing without release of lateral lid traction.

5 Per box



Muscle Hook



**BM170** Muscle Hook

5 Per box



Muscle Vectis



**BM171** Muscle Vectis

5 Per box



Scleral Depressor



**BM181** Scleral Depressor T

5 Per box

