


# SUPER VIEW™

## Tornambe Contact Lens



Superior clarity & image quality

For high magnification macular viewing

Flat vitrectomy lens offers 36° view of posterior pole

Single use for a crystal clear view every time

Economical & convenient

Product # 40413

Sold 10 per box



Insight is seeking worldwide distributors. Please contact Insight Instruments, Inc. for more information.



Insight Instruments, Inc. / 2580 SE Willoughby Blvd, Stuart, FL 34994 / Phone: 772-219-9393 / 800-255-8354  
Fax: 772-219-9342 / [www.insightinstruments.com](http://www.insightinstruments.com) / [CustomerService@insightinstruments.com](mailto:CustomerService@insightinstruments.com)

B-40413 R3 ECN-935