


Vitreoretinal Accessories


VFFBA Acar Foreign Body Forceps®

This foreign body forceps has a basket end especially designed for round objects. This end gets out from the port after entering into the eye which is made for easy entry.



VFKFB Foreign Body Forceps

This foreign body forceps has three ports which is extendible and gets out from the port after entering into the eye. The three ports help to catch the foreign body precisely and prevents disconnection once it is grasped.



VB41EA Extendible 41g Subretinal Injection Needle



VPS23 Extendible Curved Spatula 23g



